microlife®

BC50

Electric Nasal Aspirator

EN	→	1	SL	→	19
RU	→	4	SR	→	22
BG	→	7	HU	→	25
RO	→	10	HR	→	28
CZ	→	13	PL	→	31
SK	→	16			

Asia Connection Co., Ltd.
16 F-1, No.16, Jian Ba Rd.
Chung Ho District, New Taipei City
235603 Taiwan


Advena Ltd.
Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013
Malta

Distributor
Microlife AG
Espenstrasse 139
9443 Widnau / Switzerland
www.microlife.com

CE1639

IB BC 50 E-V11 1522
Revision Date: 2022-03-31

Name of Purchaser / Ф.И.О. покупателя /
 Име на купувача / Numele cumpărătorului / Imię
 i nazwisko nabywcy / Jméno kupujícího / Meno
 zákazníka / Ime in priimek kupca / Ime i prezime
 kupca / Vásárló neve / Ime i prezime kupca

Serial Number / Серийный номер / Серийен
 номер / Număr de serie / Numer serijny /
 Výrobní číslo / Výrobné číslo / Serijska številka /
 Serijski broj / Sorozatszám / Serijski broj

Date of Purchase / Дата покупки / Дата на
 закупуване / Data cumpărării / Data zakupu /
 Datum nákupu / Dátum kúpy / Datum nakupa /
 Datum kupovine / Vásárlás dátuma / Datum
 kupovine

Specialist Dealer / Специализированный дилер /
 Специалист дистрибутор / Distribuitor de spe-
 cialitate / Przedstawiciel / Specializovaný dealer /
 Špecializovaný predajca / Specializirani trgovec /
 Ovlašteni diler / Forgalmazó / Ovlašteni prodavač

- ① ON/OFF button
- ② Mucus cup
- ③ Collar
- ④ Upper cover
- ⑤ Lower cover
- ⑥ Silicone tip (size S and L)
- ⑦ Silicone rings
- ⑧ Inserting the batteries
- ⑨ Preparation and usage of this device
- ⑩ Disassembling the device
- ⑪ Washing the parts
- ⑫ Reassembling the device
- ⑬ Cleaning the mucus cup

Dear Customer,

This Microlife electric nasal aspirator provides gentle and continuous suction to remove nasal mucus easily and effectively.

If you have any questions, problems or want to order spare parts please contact your local Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the internet at www.microlife.com where you will find a wealth of invaluable information on our products.

Retain instructions in a safe place for future reference.

Stay healthy – Microlife AG!

Table of contents

1. Explanation of Symbols
2. Important Safety Instructions
3. Directions for Use
4. Inserting the batteries
5. Cleaning and Disinfecting

6. Guarantee

7. Technical Specifications

1. Explanation of Symbols

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Read the instructions carefully before using this device.

Type BF applied part

Manufacturer

Lot number

IP22

Protected against dripping water

Authorized representative
in the European Community

Operating conditions:
5 - 40 °C / 41 - 104 °F

Storage conditions:
-10 - +50°C / 14 - 122 °F

CE 1639

CE Marking of Conformity

2. Important Safety Instructions

- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.
- Prior to using the device for the first time, we recommend cleaning it as described in the section «Cleaning and Disinfecting»
- Do not use this device if you think it is damaged or notice anything unusual.
- Never open this device.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Protect it from:
 - water and moisture
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight
 - heat and cold
- The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m (according to 60601-1-2 table 5). In cases where you suspect this to be unavoidable, please verify if the device is working properly before use.
- If the device is not going to be used for a prolonged period the batteries should be removed.

Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

This device is not a toy! Keep out of reach of children and do not leave children alone during treatment.

3. Directions for Use

- Choose the suitable silicone tip ⑥ based on the size of patient's nostrils.
- Be sure to place the silicone tip ⑥ securely on the mucus cup.
- Gently insert the silicone tip into patient's nostril. Do not insert it too deeply.
- Press the ON/OFF button ① for a few seconds to collect mucus. Stop pressing the ON/OFF button when finished using the device and then the power will be stopped.

If the accumulated volume of mucus is nearly over the safe capacity of the mucus cup ② (about 4 ml), follow the section of «5. Cleaning and Disinfecting» to clean the device first.

Do not use the device if the mucus level is nearly over the opening channel (where the arrow points to on the right ⑬-a).

4. Inserting the batteries

After you have unpacked your device, first insert the batteries. Remove the battery cover by separating the upper cover ④ away from the lower cover ⑤ in a straight direction. Do not bend the covers by brute force as it can be slid open easily. Insert the batteries (2 x 1.5 V, size AA), thereby observing the indicated polarity. Check to make sure that the two silicone rings ⑦ are properly seated on the grooves. Close the covers until the "click" sound is heard.

Replace the batteries when the suction power becomes weak during use.

Always replace both batteries at the same time.

5. Cleaning and Disinfecting

- Be sure that the upper ④ and lower ⑤ covers are closed firmly before cleaning the device ⑩-a.
- Remove the silicone tip ⑩-b.
- Hold the mucus cup ② with the part between thumb and index finger ⑩-c.
- Bend the mucus cup ② downwards gently to remove it from the collar ⑩-d.
- Remove the silicone rings ⑦ from the collar ③ for cleaning after several uses of the device ⑩-e. The silicone rings do not need to be cleaned every time.

- Wash the parts thoroughly with clean water and soap (11)-a. Rinse the collar (3) under the tap (11)-b. Dry the parts before assembling them.
- Reassemble all parts (12). Properly align the positioning hole on the mucus cup (12)-b.
- Be sure that all parts fit tightly together on the device.

When the accumulated volume of mucus is over the safe capacity of the mucus cup during use of the device, if you see mucus dripping out of the air vent on the back side of the collar (13)-a, please follow instructions below:

- Be sure that the upper (4) and lower (5) covers are closed firmly (13)-b.
- Remove the silicone tip (6) and mucus cup (2).
- Have the collar (3) immersed in clean water. Press the ON/OFF button (1) for 15 seconds and the inner pump will be cleaned (13)-c.
- The water flows in and out of the inner pump in the direction as shown (13)-d.
- Remove the device from water once the cleaning is completed. Keep pressing the ON/OFF button (1) for approximately 30 seconds until the inner pump is dried.
- Dry the collar (3) and reassemble the parts (12).

6. Guarantee

This device is covered by a **3 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge. Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Battery.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website:

www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

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Electromagnetic emissions

The BC 50 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 50 should assure that it is used in such an environment.

Emissions	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The BC 50 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The BC 50 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic immunity

The BC 50 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 50 should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines	Not applicable	Not applicable
Surge IEC 61000-4-5	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5s	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	If image distortion occurs, it may be necessary to position the Model BC 50 image intensifier further from sources of power frequency magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

Electromagnetic immunity

The BC 50 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 50 should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the BC 50, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Recommended separation distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BC 50 is used exceeds the applicable RF compliance level above, the BC 50 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BC 50.